
Certificate Of Analysis

Product Information

Part Number: 35653-10
Description: OAKTON 447uS/CM COND POUCHES, 20 PACK
Lot Number: TQ2A
Expiration Date: 10/31/2017
Method Of Analysis: CONDUCTIVITY

Actual Analysis

Acceptance Range	Results	Unit Of Measure
447 US/CM +/- 5 US/CM AT 25 C	448	US/CM AT 25 C

Traceability

N.I.S.T Chemical	SRM#
POTASSIUM CHLORIDE	918b

Certified By: Thermo Fisher Scientific
Water and Lab Products
Quality Assurance

Certificate Date: 11/02/15

thermoscientific.com/water

© 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

Water and Lab Products

North America
Toll Free: 1-800-225-1480
Tel: 1-978-232-6000
info.water@thermofisher.com

Germany
Tel: (49) 6184-90-6000
info.water.uk@thermofisher.com

India
Tel: (91) 22-4157-8800
wai.asia@thermofisher.com

Japan
Tel: (81) 045-453-9175
wai.asia@thermofisher.com

China
Tel: (86) 21-68654588
wai.asia@thermofisher.com

Singapore
Tel: (65) 6778-6876
wai.asia@thermofisher.com

Australia
Tel: (613) 9757-4300
in Australia (1300) 735-295
InfoWaterAU@thermofisher.com

Thermo
S C I E N T I F I C

A Thermo Fisher Scientific Brand